



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® G Starter Kit, with P20G, P200G, P1000G microliter pipettes (metal ejector) and pipette tips

€794.00
plus VAT &
Shipping

Product Images



Description

Inexpensive in one set

Order No. 7-4423

Pipetman G Starter Kit

- Gilson Pipetman G P20G, P200G and P1000G
- Tips D200 and D1000 in Tippack Rack
- 3 single holders
- 3 plastic tip ejectors
- 1 Gilson Guide to Pipetting

Additional Information

No.	CF-0415
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4423
EAN	4058072221874
Transport temperature	Room temperature
autoclavable at 121°C	No
Number of channels	1
Functionality	mechanical
Model Microliter pipette/ Dispenser	Pipetman G P20G Pipetman G P200G Pipetman G P1000G
Type Microliter pipette (design)	Air cushion pipette

