



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson PIPETMAN® L Starter Kit, with P20L, P200L, P1000L microliter pipettes (metal ejector) and pipette tips

€940.00
plus VAT &
Shipping

Product Images

NEU

3 pipetman® L

3 SINGLE HALTER

2 TIPACK RACKS

1 PIPETTING GUIDE

...und mehr

Description

Cost effective in one kit

Gilson Pipetman® L Starter Kit:

- Gilson Pipetman® L P20L, P200L, P1000L
- Tips D200 and D1000 in Tipack Rack
- 3 single holders
- 3 plastic dispensers
- 1 User Guide Pipetman L
- 1 Two-Minute Inspection Poster
- 1 Gilson Guide to Pipetting

Additional Information

No.	CF-0406
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4443
EAN	4058072162689
Transport temperature	Room temperature
autoclavable at 121°C	No
Number of channels	1
Functionality	mechanical
Model Microliter pipette/ Dispenser	Pipetman L P20L Pipetman L P200L Pipetman L P1000L
Type Microliter pipette (design)	Air cushion pipette

