



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex TwiXS Pack E, two Acura® manual 826 XS laboratory micropipettes, with volume 0.1 – 2 µl / 2 – 20 µl

**€559.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Pipette packages Acura® manual XS 826 TwiXS Pack** for a wide range of volume combinations up to 1000 µl. The low-cost alternative to buying them individually.

Socorex Acura manual microliter pipettes have been designed with the highest customer demands in mind, ensuring convenience and reliable results for both routine and complicated pipetting processes.

The **Acura manual 826XS models** are ideally **suited for research laboratories**.

The short shaft allows optimized instrument guidance. Ideal for precise applications, such as pipetting into microtubes. The pipette line is equipped with a precision adjustment mechanism and is continuously adjustable up to a very fine volume resolution.

The innovative sealing ring concept\* stands for extremely smooth pipetting and reduced hand fatigue during workflows.

The extremely smooth volume adjustment allows simultaneous one-hand operation. The ejector is adjustable in length. Sensitive click-stops, together with the freely rotating Smartie push-button cap, ensure extreme stability of the volume setting during pipetting operations.

- Excellent ergonomics
- Smooth operation of all functions
- Unsurpassed performance and durability
- Adjustable ejection\* for the most common tips
- Integrated In-Lab calibration system Swift-Set\*
- Shock and UV light resistant, autoclavable
- CE certified IVD 98/79 EC
- Three years warranty

\* Socorex patented

## Additional Information

|                                     |   |
|-------------------------------------|---|
| No.                                 | SR-0096   |
| Manufacturer (Brand)                | Socorex   |
| old neolab article no.              | 7-7307  |
| EAN                                 | 4058072162993   |
| Transport temperature               | Room temperature  |
| autoclavable at 121°C               | No  |
| Microliter pipettes features        | IVD Compliant   |
| Number of channels                  | 1   |
| Functionality                       | mechanical  |
| Model Microliter pipette/ Dispenser | Socorex Acura manual  |
| Type Microliter pipette (design)    | Air cushion pipette   |
| Scope of delivery                   | Two Acura® manual XS pipettes, shelf pipette holder, Qualitips® tip sample, control certificate, operating instructions |

