



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand Haematocrit sealing kit, 24 storage fields for 2 capillaries each, CE-IVD

€42.20
plus VAT &
Shipping

Product Images



Description

Non-drying putty on a plastic plate. For rapid sealing of micro-hematocrit capillaries. 2 capillaries each (24 in total) can be placed vertically in the numbered storage fields on the long sides of the plate. CE-marked in accordance with IVD Directive 98/79 EC. Note:

Hold the capillary, which is approximately 2/3 full, close to the empty end and, to close it, insert the empty end into the putty, turning it slightly.

Additional Information

No.	BN-2084
Manufacturer (Brand)	Brand
EAN	4058072175740
Gross weight	0.1824
HSNumber	32141010
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.16120 kg
Country of origin	Slovakia
Packing width	0.08900 m
Packing height	0.03400 m
Packing depth	0.12100 m
Packaging volume	366.14600000 cc
Customs tariff number	32141010

