



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hirschmann® disposable capillary pipettes with ring mark at 44.7 µl NA-hep, 250 pcs/pack

**Price on
request
plus VAT &
Shipping**



Description

First class quality capillaries, made of high precision capillary tube.
Smooth, clean surfaces with clear color code according to DIN, ring mark and volume.
Lower end cut cleanly and exactly at right angles so that the liquid column breaks off in a defined manner.

Order No. 1-7306

Sodium-heparinized (225 i.E. sodium)

Accessories: Micropipetter (Order No. 1-7315)

Additional Information

No.	NH-1653
Manufacturer (Brand)	Hirschmann
old neoLab article no.	1-7306
EAN	4058072336776
Transport temperature	Room temperature
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 1
Volume	44.7 µl
Length	130 mm
Type pipette	Capillary pipettes

