



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand microhematocrit capillaries, surface heparinized, color code red, length 75mm, outer- $\varnothing$  1,55 mm**

**€30.80** As low  
as  
**€29.20**  
plus VAT &  
Shipping

## Product Images



## Description

---

DIN ISO 12772 and BS 4316-68. With color code. For use in hematocrit centrifuges. Heparinized: Entire inner surface with sodium heparinate, with red color code, for capillary blood collection. Non-heparinized: With blue color code, for determination from heparinized venous blood. CE-marked according to IVD Directive 98/79 EC. DIN ISO 12772 and BS 4316-68. With color code. For use in hematocrit centrifuges. Heparinized: Entire inner surface with sodium heparinate, with red color code, for capillary blood collection. Non-heparinized: With blue color code, for determination from heparinized venous blood. CE-marked according to IVD Directive 98/79 EC.

## Additional Information

---

No.	BN-2082
Manufacturer (Brand)	Brand
EAN	4058072177249
Gross weight	0.3131
HSNumber	70172000
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.16000 kg
diameter	1,55 mm
Country of origin	Denmark
Packing width	0.08200 m
Packing height	0.04700 m
Packing depth	0.11600 m
Packaging volume	447.06400000 cc
Customs tariff number	70172000
Length	75 mm
Color code	red

