



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand microhematocrit capillaries, surface not heparinized, color code blue, length 75mm, outer- \varnothing 1,55 mm

€29.35 As low
as
€27.40
plus VAT &
Shipping

Product Images



Description

DIN ISO 12772 and BS 4316-68. With color code. For use in hematocrit centrifuges. Heparinized: Entire inner surface with sodium heparinate, with red color code, for capillary blood collection. Non-heparinized: With blue color code, for determination from heparinized venous blood. CE-marked according to IVD Directive 98/79 EC.

Additional Information

No.	BN-2083
Manufacturer (Brand)	Brand
VGKL number	144974932
EAN	4058072181376
Gross weight	0.3262
HSNumber	70179000
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.30600 kg
diameter	1,55 mm
Country of origin	Denmark
Packing width	0.08300 m
Packing height	0.05000 m
Packing depth	0.11600 m
Packaging volume	481.40000000 ccm
Customs tariff number	70179000
Length	75 mm
Color code	blue

