



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Hirschmann® graduated pipette 25.0:0.1 ml
wide aperture, clear glass, graduated, split to
tip**

€27.42 As low as
€22.15
**plus VAT &
Shipping**

Product Images



Description

Graduated pipette, Class B, AR® glass, permanently graduated in brown with silver diffusion dye, wide aperture. Serology, line graduation, with extended tip opening. Zero point at top, divided to tip. From 5 ml with cotton stopper end as standard.

Additional Information

No.	NH-0017
Manufacturer (Brand)	Hirschmann
old neoLab article no.	E-0414
EAN	4058072336011
Transport temperature	Room temperature
Color	transparent (clear)
Material	Borosilicate glass
Volume	25 ml
Length	450 mm
Volumetric instruments features	with cotton plug end
Color graduation	brown
Class	B
Position zero point	Top
Division	0.1 ml
Type pipette	Serological pipette with wide tip opening
adjusted to	"ex"

