



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Poulsen & Graf Glass Pasteur pipettes with  
cotton stopper, non-sterile, 50 mm lg., 250  
pcs./pack**

**€27.59  
plus VAT &  
Shipping**

## Product Images



## Description

---

With finely extended tip.  
ISO 7712, all dimensions are approximate.

### **Order No. E-1098**

with cotton plug

## Additional Information

---

No.	PL-1989
Manufacturer (Brand)	Poulsen & Graf
old neoLab article no.	E-1098
EAN	4058072338053
Transport temperature	Room temperature
diameter	7 mm
Color	transparent (crystal clear)
Material	Soda-lime glass (normal glass) 1
Volume	2 ml
sterile	No
Length	150 mm
Volumetric instruments features	with cotton plug end
Drops per ml	50
Drop size	20 µl
Type pipette	Pasteur pipettes

