



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Poulten & Graf Glass Pasteur pipettes with
cotton stopper, non-sterile, 230 mm lg., 250
pcs./pack**

**€30.44
plus VAT &
Shipping**

Product Images



Description

With finely extended tip.
ISO 7712, all dimensions are approximate.

Order No. E-1099
with cotton stopper

Additional Information

No.	PL-1994
Manufacturer (Brand)	Poulsen & Graf
old neoLab article no.	E-1099
EAN	4058072338077
Transport temperature	Room temperature
diameter	7 mm
Color	transparent (crystal clear)
Material	Soda-lime glass (normal glass) 1
Volume	2 ml
sterile	No
Length	230 mm
Volumetric instruments features	with cotton plug end
Drops per ml	50
Drop size	20 µl
Type pipette	Pasteur pipettes

