



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Poulten & Graf Pasteur pipette, 230 mm, glass with cotton stopper, sterile

€11.90
plus VAT &
Shipping

Product Images



Description

The sterile neoLab Pasteur pipettes are made of soda-lime glass and comply with the ISO 7712 standard. With cotton stopper.

Additional Information

No.	PL-1996
Manufacturer (Brand)	Poulsen & Graf
old neoLab article no.	1-0206
EAN	4058072338084
Transport temperature	Room temperature
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 1
sterile	Yes
Length	230 mm
Volumetric instruments features	with cotton plug end
Type pipette	Pasteur pipettes
standardized according to ISO standard	ISO 7712

