



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Pasteur pipettes 7.3 ml, ungraduated,
300 mm long, 100 pcs./pack**

**€41.30
plus VAT &
Shipping**

Product Images



Description

High quality safety pasteur pipettes

- made of one piece
- without seam
- with integrated suction ball, which can be easily pressed and still develops full suction power.

Fatigue-free work even with longer pipetting series.

CE-marked according to IVDD 98/79/EEC (Declaration of Conformity available on request).

Additional Information

No.	1-6163
Manufacturer (Brand)	neoLab
EAN	4058072036966
diameter	8,7 mm
Color	transparent (milky)
Material	Low-density polyethylene (LDPE)
Volume	7.3 ml
sterile	No
Length	300 mm
Volumetric instruments features	CE marked
Drops per ml	20
Drop size	50 µl
Type pipette	Pasteur pipettes

