



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand Pasteur pipette, PE-LD, 3 ml, grad.
0.5/1.0/1.5/2.0/2.5/3.0 ml, \varnothing tip 3 mm, suction
volume max. 6.2 ml**

€37.98 As low
as
€34.28
**plus VAT &
Shipping**

Product Images



Description

PE-LD. Very good reproducibility of drop count per milliliter, therefore ideal for dispensing aliquots of liquid. Pasteur pipettes can be deep-frozen when full or converted into a closed vessel by heat sealing if required. The integrated suction bellows can be easily compressed. This keeps fingers fatigue-free even with frequent pipetting. Stable against gas or gamma ray sterilization procedures.

Additional Information

No.	BN-2078
Manufacturer (Brand)	Brand
VGKL number	284947760
EAN	4058072176730
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Country of origin	Italy
Customs tariff number	39269097
Division	0.5 ml
Drops per ml	19

