



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand Pasteur pipette, PE-LD, ungraduated, \emptyset tip
2.5 mm, suction volume max. 3.3 ml**

€47.13 As low as
€42.53
**plus VAT &
Shipping**

Product Images



Description

PE-LD. Very good reproducibility of drop count per milliliter, therefore ideal for dispensing aliquots of liquid. Pasteur pipettes can be deep-frozen when full or converted into a closed vessel by heat sealing if required. The integrated suction bellows can be easily compressed. This keeps fingers fatigue-free even with frequent pipetting. Stable against gas or gamma ray sterilization procedures.

Additional Information

No.	BN-2076
Manufacturer (Brand)	Brand
VGKL number	104947750
EAN	4058072174125
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Country of origin	Czech Republic
Customs tariff number	39269097
Drops per ml	27

