



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Moonlab® Serological Pipettes 1 ml, individually sterile packed, 500 pcs./pack

Special Price


**€55.30** was

**€79.00**

**plus VAT &  
Shipping**

## Product Images

---



## Description

---

Moonlab® brand crystal clear serological pipettes are manufactured from high quality polystyrene in a clean room environment. In addition, the serological pipettes are sterilized and individually packaged in easy-to-open foil/paper peel bags. This prevents contamination of your samples with unwanted substances (RNases, DNases, pyrogens). In the mouthpiece, cotton swabs provide protection against unwanted contamination of the pipetting aid.

- Color coding facilitates identification of the correct pipette
- Larger pipetting volume due to negative scaling
- Universal mouthpiece for precise insertion into pipetting aids

## Additional Information

---

No.	4-0130
Manufacturer (Brand)	Moonlab
EAN	4058072154684
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	1 ml
sterile	Yes
Length	300 mm
Volumetric instruments features	pyrogen-free with cotton wool/filter plug with color code graduated for 1x use
Color code	yellow
Color graduation	black
Type serological pipette	Standard

