



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Moonlab® Serological Pipettes 10 ml, individually sterile packed, 200 pcs./pack


Special Price

€44.10 was

€63.00

**plus VAT &
Shipping**

Product Images



Description

Moonlab® brand crystal clear serological pipettes are manufactured from high quality polystyrene in a clean room environment. In addition, the serological pipettes are sterilized and individually packaged in easy-to-open foil/paper peel bags. This prevents contamination of your samples with unwanted substances (RNases, DNases, pyrogens). In the mouthpiece, cotton swabs provide protection against unwanted contamination of the pipetting aid.

- Color coding facilitates identification of the correct pipette
- Larger pipetting volume due to negative scaling
- Universal mouthpiece for precise insertion into pipetting aids

Additional Information

No.	4-0133
Manufacturer (Brand)	Moonlab
EAN	4058072154196
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	10 ml
sterile	Yes
Length	340 mm
Volumetric instruments features	pyrogen-free with cotton wool/filter plug with color code graduated for 1xuse
Color code	orange
Color graduation	black
Type serological pipette	Standard

