




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Serological Pipettes 5.0 ml, sterilized


Special Price

**€41.30** was

**€59.00**

**plus VAT &  
Shipping**

## Product Images



## Description

---

- Serological pipettes made of crystal clear polystyrene (PS).
- RNase-, DNase-, Pyrogen-free - non-toxic, TSE/BSE free
- imprinted, black, permanent graduation
- with color-coded ring for easy identification, colors according to standard code
- individually sterile packed - packaging: foil/paper
- additional capacity due to negative graduation
- accuracy +/- 2-3% on full volume

## Additional Information

---

No.	DK-0934
Manufacturer (Brand)	neoLab
EAN	4066292687820
Transport temperature	Room temperature
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	5 ml
Length	348.1 mm
Volumetric instruments features	sterile pyrogen-free with color code graduated
Color code	blue
Color graduation	black
Position zero point	Top
Division	0.1 ml
Type pipette	Serological pipette with wide tip opening
Type serological pipette	Standard
adjusted to	"in"

