




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Serological Pipettes 10 ml, sterilized

Special Price

€41.30 was

€59.00

**plus VAT &
Shipping**

Product Images

Description

- Serological pipettes made of crystal clear polystyrene (PS).
- RNase-, DNase-, Pyrogen-free - non-toxic, TSE/BSE free
- imprinted, black, permanent graduation
- with color-coded ring for easy identification, colors according to standard code
- individually sterile packed - packaging: foil/paper
- additional capacity due to negative graduation
- accuracy +/- 2-3% on full volume

Additional Information

No.	DK-0935
Manufacturer (Brand)	neoLab
EAN	4066292552753
Transport temperature	Room temperature
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	10 ml
Length	346.8 mm
Volumetric instruments features	sterile pyrogen-free with color code graduated
Color code	orange
Color graduation	black
Position zero point	Top
Division	0.1 ml
Type pipette	Serological pipette with wide tip opening
Type serological pipette	Standard
adjusted to	"in"

