




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Serological Pipettes 50 ml, sterilized

€121.00
plus VAT &
Shipping

Product Images

Description

- Serological pipettes made of crystal clear polystyrene (PS).
- RNase-, DNase-, Pyrogen-free - non-toxic, TSE/BSE free
- imprinted, black, permanent graduation
- with color-coded ring for easy identification, colors according to standard code
- individually sterile packed - packaging: foil/paper
- additional capacity due to negative graduation
- accuracy +/- 2-3% on full volume

Additional Information

No.	DK-0937
Manufacturer (Brand)	neoLab
EAN	4058072997144
Transport temperature	Room temperature
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	50 ml
Length	368.2 mm
Volumetric instruments features	sterile pyrogen-free with color code graduated
Color code	black
Color graduation	black
Position zero point	Top
Division	1 ml
Type pipette	Serological pipette with wide tip opening
Type serological pipette	Standard
adjusted to	"in"

