



neoLabLine dip pipette 10 ml

Product Images



Description

For sampling from large containers, for quick and easy removal of firmly defined volumes.
Immerse pipettes and allow the hollow body to fill. The sample is emptied by tilting.

Additional Information

No.	SE-0961
Manufacturer (Brand)	neoLabLine
old neoLab article no.	2-5701
EAN	4058072004798
diameter	17 mm
Color	transparent (klar)
Material	Kalk-Natron-Glas (Normalglas) 1
Volume	10 ml
Length	270 mm
Type liquid collector	Tauchpipette
Type sampler	Flüssigkeitssammler
for medium	Flüssigkeiten

