



neoLabLine dip pipette 25 ml

Product Images



Description

For sampling from large containers, for quick and easy removal of firmly defined volumes.
Immerse pipettes and allow the hollow body to fill. The sample is emptied by tilting.

Additional Information

No.	SE-0977
Manufacturer (Brand)	neoLabLine
old neoLab article no.	2-5703
EAN	4066292683624
diameter	28 mm
Color	transparent (klar)
Material	Kalk-Natron-Glas (Normalglas) 1
Volume	25 ml
Length	305 mm
Type liquid collector	Tauchpipette
Type sampler	Flüssigkeitssammler
for medium	Flüssigkeiten

