



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand PD-Tips II single packed sterile BIO-CERT, 2.5 ml

**€167.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The **dispenser tips PD-Tips II** for efficient, step-by-step dispensing are the consistent further development of the proven BRAND PD-Tips.

Thanks to the patented coding, BRAND <sup>HandyStep®</sup> touch and <sup>HandyStep®</sup> electronic automatically recognize the volume of the tips.

The direct displacement principle enables precise work even with viscous media or those with high vapor pressure. With the tightly sealing plunger, the PD-Tips II effectively protect your dispenser from aerosols. This protects against contamination and ensures a long service life of the equipment.

- Volumes from 0.1 ml to 50 ml
- Compatible with automatic size recognition: <sup>HandyStep®</sup> touch, <sup>HandyStep®</sup> electronic
- Compatible without automatic size detection: <sup>HandyStep®</sup> S, <sup>HandyStep®</sup>
- Available in <sup>BIO-CERT®</sup> quality (thus free of DNA, RNase, endotoxins and ATP)
- Conform to the requirements of ISO 8655
- Supplied with batch certificate
- Manufactured from high-quality starting materials (cylinder: PP / piston: PE-HD, 0.1 ml: LCP)
- Adapters for 25 ml and 50 ml sizes are supplied
- PD-Tips II are not autoclavable, the adapters are autoclavable

## Additional Information

No.	BN-2636
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Brand
EAN	4058072217594
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	2500 µl
sterile	Yes
Pipette tips properties	graduated with certification DNA-free RNase-free Endotoxin-free ATP-free
Fluid behavior	hydrophobic
Packing units tips	100 x 1
for model pipette	HandyStep S HandyStep electronic Sunlab Dispenser
Type certificate	Batch certificate from the manufacturer
Packaging type pipette tips (BRAND)	individually packed
for pipettes / dispensers of the manufacturers	Brand Sunlab

