



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand PD-Tips II, loose, non-sterile, 1.25 ml

€78.25
plus VAT &
Shipping

Product Images



Description

The **dispenser tips PD-Tips II** for efficient, step-by-step dispensing are the consistent further development of the proven BRAND PD-Tips.

Thanks to the patented coding, BRAND ^{HandyStep®} touch and ^{HandyStep®} electronic automatically recognize the volume of the tips.

The direct displacement principle enables precise work even with viscous media or those with high vapor pressure. With the tightly sealing plunger, the PD-Tips II effectively protect your dispenser from aerosols. This protects against contamination and ensures a long service life of the equipment.

- Volumes from 0.1 ml to 50 ml
- Compatible with automatic size recognition: ^{HandyStep®} touch, ^{HandyStep®} electronic
- Compatible without automatic size detection: ^{HandyStep®} S, ^{HandyStep®}
- Also available in ^{BIO-CERT®} quality (thus free of DNA, RNase, endotoxins and ATP)
- Conform to the requirements of ISO 8655
- Supplied with batch certificate
- Manufactured from high-quality starting materials (cylinder: PP / piston: PE-HD, 0.1 ml: LCP)
- Adapters for 25 ml and 50 ml sizes are supplied
- PD-Tips II are not autoclavable, the adapters are autoclavable

Additional Information

No.	BN-2624
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Brand
EAN	4058072171407
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	1250 µl
sterile	No
Pipette tips properties	graduated with certification
Fluid behavior	hydrophobic
Packing units tips	1 x 100
for model pipette	HandyStep S HandyStep electronic Sunlab Dispenser
Type certificate	Batch certificate from the manufacturer
Packaging type pipette tips (BRAND)	Bag
for pipettes / dispensers of the manufacturers	Brand Sunlab

