



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Combitips® advanced, Biopur®, 0.2 mL, light blue, colorless tips, 100 pieces, individually blistered

**€195.00
plus VAT &
Shipping**

Product Images



Description

The ideal system component for all multipipette pipettes from Eppendorf

With the introduction of the Combitips® advanced/Multipipette system in 1978, Eppendorf revolutionized the industry and laboratory procedures and has been the market leader in terms of product innovation for almost 40 years. The Combitips® advanced have been completely redesigned and optimized all around to even better meet the demands of modern laboratory routine. Combitips® advanced also work according to the direct displacement principle. This means that the correct desired volume is always dispensed, regardless of the density and flow properties of the liquid (e.g. higher vapor pressure or increased viscosity). Working with radioactive or toxic substances is also safer thanks to the hermetically sealed piston, as aerosol contamination is prevented.

Additional Information

| | |
|-----------------------|--|
| No. | VB-0256 |
| Article - USP | Combitips® advanced are the ideal system component for all Multipette pipettes from Eppendorf Available in nine volume sizes (0.1 mL - 50 mL), which can be used to dispense a wide range of volumes, depending on the Combitips® advanced used |
| Typ Pipettenspitze | Direktverdrängerspitzen |
| Manufacturer (Brand) | Eppendorf |
| VGKL number | 106989626 |
| EAN | 4058072172589 |
| Gross weight | 0.63 |
| HSNumber | 84799070 |
| Storage temperature | Room temperature |
| Transport temperature | Room temperature |
| autoclavable at 121°C | No |
| Net weight | 0.56700 kg |
| Volume MAX | 200 µl |
| Volume MIN | 40 µl |
| Country of origin | Germany |
| sterile | Yes |
| Packing width | 0.42000 m |
| Packing height | 0.16000 m |
| Packing depth | 0.17000 m |
| Packaging volume | 12000.000000 ccm |
| Customs tariff number | 84799070 |
| Degree of purity | Biopur |

