



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Combitips® advanced, Forensic DNA Grade, 5.0 mL, blue, colorless tips, 100 pieces, individually blistered

€150.00
plus VAT &
Shipping

Product Images



Description

The ideal system component for all multipipette pipettes from Eppendorf

With the introduction of the Combitips® advanced/Multipipette system in 1978, Eppendorf revolutionized the industry and laboratory procedures and has been the market leader in terms of product innovation for almost 40 years. The Combitips® advanced have been completely redesigned and optimized all around to even better meet the demands of modern laboratory routine. Combitips® advanced also work according to the direct displacement principle. This means that the correct desired volume is always dispensed, regardless of the density and flow properties of the liquid (e.g. higher vapor pressure or increased viscosity). Working with radioactive or toxic substances is also safer thanks to the hermetically sealed piston, as aerosol contamination is prevented.

Additional Information

No.	VB-0280
Article - USP	Combitips® advanced are the ideal system component for all Multipette pipettes from Eppendorf Available in nine volume sizes (0.1 mL - 50 mL), which can be used to dispense a wide range of volumes, depending on the Combitips® advanced used
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072187323
Gross weight	0.93
HSNumber	84799070
Storage temperature	Room temperature
Transport temperature	Room temperature
autoclavable at 121°C	No
Net weight	0.83700 kg
Volume MAX	5 ml
Volume MIN	1000 µl
Country of origin	Germany
sterile	Yes
Packing width	0.46000 m
Packing height	0.16000 m
Packing depth	0.19000 m
Packaging volume	14000.000000 ccm
Customs tariff number	84799070
Degree of purity	Forensic DNA Grade

