



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Combitips advanced, 0.1 ml, white, 100 pcs./pack

€107.00
plus VAT &
Shipping

Product Images



Description

Combitips^{advanced®} for use in all Eppendorf^{Multipettes®}

- individual color code
- extended tips provide more flexibility
- new shape allows more ergonomic insertion of CBT advanced
- optimized printing
- adapter for 25/50 ml included (Order No. E-6407/08)
- new, resealable packaging / dispenser box. (4 bags à 25 pieces)
- made of pure PP, without lubricants, biocides and plasticizers
- work according to the direct displacement principle, and therefore suitable for working with liquid with high viscosity

Additional Information

No.	VB-0245
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6400
EAN	4058072405328
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	100 µl
Volume	0.1 ml
sterile	No
Pipette tips properties	with certification graduated
Fluid behavior	hydrophobic
Manufacturer	Eppendorf
Type certificate	Batch certificate from the manufacturer
Packaging type pipette tips (general)	Cardboard
Packing units tips	1 x 100
for pipettes / dispensers of the manufacturers	Eppendorf

