



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Combitips advanced, 5.0 ml, blue, 100 pcs./pack

€107.00
plus VAT &
Shipping

Product Images



Description

Combitips advanced® for use in all Eppendorf Multipettes®

- individual color code
- extended tips provide more flexibility
- new shape enables more ergonomic insertion of CBT advanced
- optimized printing
- adapter for 25/50 ml included (Order No. E-6407/08)
- new, resealable packaging / dispenser box. (4 bags à 25 pieces)
- made of pure PP, without lubricants, biocides and plasticizers
- work according to the direct displacement principle, and therefore suitable for working with liquid with high viscosity

Additional Information

No.	VB-0250
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6405
EAN	4058072405373
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	5 ml
Volume	5 ml
sterile	No
Pipette tips properties	with certification graduated
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Cardboard
Packing units tips	1 x 100
Manufacturer	Eppendorf
Type certificate	Batch certificate from the manufacturer
for pipettes / dispensers of the manufacturers	Eppendorf

