



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Combitips advanced, 25 ml, red, 100 pcs./pack

€126.00
plus VAT &
Shipping

Product Images



Description

Combitips^{advanced®} for use in all Eppendorf^{Multipettes®}

- individual color code
- extended tips provide more flexibility
- new shape enables more ergonomic insertion of CBT advanced
- optimized printing
- adapter for 25/50 ml included (Order No. E-6407/08)
- new, resealable packaging / dispenser box. (4 bags à 25 pieces)
- made of pure PP, without lubricants, biocides and plasticizers
- work according to the direct displacement principle, and therefore suitable for working with liquid with high viscosity

Order No. E-6407: 25 ml/red

Additional Information

No.	VB-0252
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6407
EAN	4058072405397
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	25 ml
Volume	25 ml
sterile	No
Pipette tips properties	with certification graduated
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Cardboard
Packing units tips	1 x 100
Spare part	VB-0264
Manufacturer	Eppendorf
Type certificate	Batch certificate from the manufacturer
for pipettes / dispensers of the manufacturers	Eppendorf

