



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Combitips advanced, 50 ml, light gray, 100 pcs./pack

**€139.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

### **Combitips<sup>advanced®</sup> for use in all Eppendorf<sup>Multipettes®</sup>**

- individual color code
- extended tips provide more flexibility
- new shape enables more ergonomic insertion of CBT advanced
- optimized printing
- adapter for 25/50 ml included (Order No. E-6407/08)
- new, resealable packaging / dispenser box. ( 4 bags à 25 pieces)
- made of pure PP, without lubricants, biocides and plasticizers
- work according to the direct displacement principle, and therefore suitable for working with liquid with high viscosity

## Additional Information

No.	VB-0253
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6408
EAN	4058072405403
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	50 ml
Volume	50 ml
sterile	No
Pipette tips properties	with certification graduated
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Cardboard
Packing units tips	1 x 100
Spare part	VB-0265
Manufacturer	Eppendorf
Type certificate	Batch certificate from the manufacturer
for pipettes / dispensers of the manufacturers	Eppendorf

