



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Combitips, advanced 25.0ml (color code: red) biopur

€230.00
plus VAT &
Shipping

Product Images



Description

Combitips® advanced, Biopur®, 25 mL, red, colorless tips, 100 pieces (4 bags × 25 pieces), individually blistered, each box complete with one adapter – Combitips® advanced are the ideal system component for all Multipipette pipettes from Eppendorf – Available in nine volume sizes (0.1 mL – 50 mL), with which a wide variety of volumes can be dispensed, depending on the Combitips® advanced used. The smallest step size is 100 nanoliters – Extended Combitips® advanced for 2.5 mL, 5.0 mL and 10 mL volumes extend to the bottom of almost all common laboratory tubes – Individual color coding facilitates quick identification of the desired Combitips® advanced – Unique funnel geometry ensures comfortable handling and prevents damage to gloves – Combitips® advanced are available in different purity grades and thus suitable for any type of application – The right purity grade for every application: Eppendorf Quality™, PCR Clean, Forensic DNA Grade or Biopur®

Additional Information

No.	VB-0262
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072172060
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	25 ml
sterile	No
Pipette tips properties	with certification graduated
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Cardboard
Packing units tips	1 x 100
Type certificate	Batch certificate from the manufacturer
for pipettes / dispensers of the manufacturers	Eppendorf

