



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Combitips, advanced 1.0ml (color code: yellow) PCR clean

**€134.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Combitips® advanced, PCR clean, 1.0 mL, yellow, colorless tips, 100 pieces (4 resealable bags × 25 pieces) – Combitips® advanced are the ideal system component for all Multipette pipettes from Eppendorf – Available in nine volume sizes (0.1 mL – 50 mL), with which a wide variety of volumes can be dispensed, depending on the Combitips® advanced used. The smallest step size is 100 nanoliters – Extended Combitips® advanced for 2.5 mL, 5.0 mL and 10 mL volumes extend to the bottom of almost all common laboratory tubes – Individual color coding facilitates quick identification of the desired Combitips® advanced – Unique funnel geometry ensures comfortable handling and prevents damage to gloves – Combitips® advanced are available in different purity grades and thus suitable for any type of application – The right purity grade for every application: Eppendorf Quality™, PCR Clean, Forensic DNA Grade or Biopur®

## Additional Information

---

No.	VB-0272
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072170912
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	1000 µl
sterile	No
Pipette tips properties	with certification graduated
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Cardboard
Packing units tips	1 x 100
Type certificate	Batch certificate from the manufacturer
for pipettes / dispensers of the manufacturers	Eppendorf

