



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Finntip Stepper, 1.25 ml, sterile, 50 pcs./box

€133.00
plus VAT &
Shipping

Product Images



Description

Pipette tips for Finntip™ stepper

The tips are individually sterile packed.

The accuracy is ± 0.5 to ± 1.0 % depending on the tip and volume.

Additional Information

No.	PO-0312
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1467
EAN	4058072338374
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	1250 µl
Volume	1.25 ml
sterile	Yes
Pipette tips properties	graduated
Fluid behavior	hydrophobic
Manufacturer	Thermo Scientific
Packaging type pipette tips (general)	Carton. Tips individually packed
Packing units tips	1 x 50
for pipettes / dispensers of the manufacturers	Thermo Scientific

