



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Finntip Stepper, 5.0 ml, sterile, 25 pcs./box

€66.80
plus VAT &
Shipping

Product Images



Description

Pipette tips for Finntip™ stepper

The tips are individually sterile packed.
The accuracy is ± 0.5 to ± 1.0 % depending on the tip and volume.

Additional Information

No.	PO-0316
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1469
EAN	4058072338398
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	5 ml
Volume	5 ml
sterile	Yes
Pipette tips properties	graduated
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Carton. Tips individually packed
Packing units tips	1 x 25
Manufacturer	Thermo Scientific
for pipettes / dispensers of the manufacturers	Thermo Scientific

