



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Thermo Scientific Finntip Stepper, 12.5 ml, sterile, 25 pcs./box

**€66.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Pipette tips for Finntip™ stepper

The tips are individually sterile packed.

The accuracy is  $\pm 0.5$  to  $\pm 1.0$  % depending on the tip and volume.

## Additional Information

---

No.	PO-0318
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	3-1470
EAN	4058072338404
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	12,5 ml
Volume	12.5 ml
sterile	Yes
Pipette tips properties	graduated
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Carton. Tips individually packed
Packing units tips	1 x 25
Manufacturer	Thermo Scientific
for pipettes / dispensers of the manufacturers	Thermo Scientific

