



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## **Gilson DistriTips for Distriman, Mini 1250 µl, 20-100 µl, 50 pcs/pack**

**€54.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

**Distritips® for use with the DISTRIMAN®**

Loosely packed standard tips, non-sterile.

- only 3 tips cover the entire volume range (aliquots from 1 µl to 1.25 ml)
- specifications exceed ISO 8655 guidelines
- the Distritips® work according to the direct displacement principle (suitable for viscous liquids)

## Additional Information

No.	CF-0358
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4383
EAN	4058072221430
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	1250 µl
Volume	1.25 ml
sterile	No
Pipette tips properties	with certification graduated
Designation tip (Gilson)	Mini
Fluid behavior	hydrophobic
Manufacturer	Gilson
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Tips and plungers in carton
Packaging type pipette tips (general)	Cardboard
Packing units tips	1 x 50
for pipettes / dispensers of the manufacturers	Gilson

