



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson DistriTips for Distriman, Maxi 12.5 ml, 200-1000 µl, 50 pcs/pack

€54.00
plus VAT &
Shipping

Product Images



Description

Distritips® for use with the DISTRIMAN®

Loosely packed standard tips, non-sterile.

- only 3 tips cover the entire volume range (aliquots from 1 µl to 1.25 ml)
- specifications exceed ISO 8655 guidelines
- the Distritips® work according to the direct displacement principle (suitable for viscous liquids)

Additional Information

No.	CF-0360
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4385
EAN	4058072221447
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	12,5 ml
Volume	12.5 ml
sterile	No
Pipette tips properties	with certification graduated
Fluid behavior	hydrophobic
Designation tip (Gilson)	Maxi
Packaging type pipette tips (general)	Cardboard
Packing units tips	1 x 50
Manufacturer	Gilson
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (Gilson)	Tips and plungers in carton
for pipettes / dispensers of the manufacturers	Gilson

