



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Gilson DistriTips for Distriman, Maxi ST 12.5 ml, 200-1000 µl, 50 pcs/pack

€106.00
plus VAT &
Shipping

Product Images



Description

Distritips® for use with the DISTRIMAN®

The Distritips® are sterilized and individually packaged. They are free of RNase, DNase, DNA, RNA and proteases.

- only 3 tips cover the entire volume range (aliquots from 1 µl to 1.25 ml)
- specifications exceed ISO 8655 guidelines
- the Distritips® work according to the direct displacement principle (suitable for viscous liquids)

Additional Information

No.	CF-0364
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4386
EAN	4058072221478
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	12,5 ml
Volume	12.5 ml
sterile	Yes
Pipette tips properties	pyrogen-free with certification graduated RNase-free RNA-free DNase-free DNA-free ATP-free
Designation tip (Gilson)	Maxi
Fluid behavior	hydrophobic
Manufacturer	Gilson
Type certificate	Sterilization certificate from the manufacturer
Packaging type pipette tips (Gilson)	Tips and plungers in carton
Packaging type pipette tips (general)	Carton. Tips individually packed
Packing units tips	1 x 50
for pipettes / dispensers of the manufacturers	Gilson

