



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Ecostep™ syringes for Stepper 411 blue 3.75 ml, range 50-500 µl, PU: 100 pcs/pack

€112.94
plus VAT &
Shipping

Product Images



Description

Ecostep™ syringes for Stepper™ 411

Three sizes cover a volume range from 10 to 5000 µl.

Color coding on syringes and selection button prevent any setting errors.

Graduation helps to control syringe contents.

PE/PP materials guarantee excellent chemical resistance.

Not sterilized; packed in bags

Additional Information

No.	SR-0152
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Socorex
old neolab article no.	7-7419
EAN	4058072163068
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (blue)
Volume MAX	3,75 ml
Volume	3.75 ml
sterile	No
Pipette tips properties	graduated
Manufacturer	Socorex
Packaging type pipette tips (general)	Cardboard
Packing units tips	1 x 100
for pipettes / dispensers of the manufacturers	Socorex

