



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Ecostep bioproof™ for Stepper 411 red 37.5 ml sterile, single, 50 p

€174.93
plus VAT &
Shipping

Product Images



Description

Ecostep™ syringes for Stepper™ 411

Three sizes cover a volume range from 10 to 5000 µl.

Color coding on syringes and selection button prevent any setting errors.

Graduation helps to control syringe contents.

PE/PP materials guarantee excellent chemical resistance.

bioproof: individually packaged, sterilized syringes with high biological purity. Each production series is tested and certified free of DNase, RNase, pyrogens (endotoxins) and human DNA by independent testing laboratories.

Additional Information

No.	SR-0157
Typ Pipettenspitze	Direktverdrängerspitzen
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7424
EAN	4058072163105
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (red)
Volume MAX	37,5 ml
Volume	37.5 ml
sterile	Yes
Pipette tips properties	pyrogen-free graduated RNase-free RNA-free DNase-free DNA-free
Manufacturer	Socorex
Type certificate	Sterilization certificate from the manufacturer
Packaging type pipette tips (general)	Carton. Tips individually packed
Packing units tips	1 x 50
for pipettes / dispensers of the manufacturers	Socorex

