



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



AHN Filter tips myTipLow Retention 100-1000 µl, clear, sterile

€93.44 As low

as

€71.88

**plus VAT &
Shipping**

Product Images



Description

AHN myTip® Low Retention Filter Tips

Thanks to the combination of the extremely hydrophobic inner surface and the micropore filter technology, the ultimate performance of the AHN myTip® Low Retention Filter Tips is guaranteed.

Precise and comfortable pipetting with maximum sample recovery has never been easier than with AHN myTip® Low Retention Filter Tips. With AHN myTip® Low Retention Filter Tips you can precisely dispense even the smallest sample volumes, ensure the reproducibility of your results and increase your working accuracy at the same time. Viscous liquids can be processed more easily with AHN myTip® Low Binding Filter Tips.

Create an effective barrier to harmful aerosol transfer while processing your samples with AHN myTip® Low Retention filter tips. Protect the inside of your pipette and reduce repair and service costs. AHN myTip® Low Retention filter tips reduce the risk of cross-contamination and ensure stable and reliable results.

AHN myTip® Low Binding Filter Tips are made of specially designed, hydrophobic, pure polyethylene that is free of chemical additives and eliminates the risk of cross-contamination should the liquid come into contact with the filter.

AHN myTip® Low Retention Filter Tips can be universally used with common brands of pipettes.

Additional Information

No.	AN-0346
Typ Pipettenspitze	Filterspitzen low retention
Manufacturer (Brand)	AHN
old neoLab article no.	7-3125
EAN	4058072083816
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume MAX	1000 µl
Volume MIN	100 µl
sterile	Yes
Pipette tips properties	with filter with box
Color code	blue
Fluid behavior	hydrophobic
Model pipette tip	AHN Low Retention Filter Tips
Packaging type pipette tips (AHN)	Rack

