



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand filter tips, PP, 5-200 µl TipBox, BIO-CERT (sterile) Packing unit: 10 x 96 pcs

Special Price

**€68.78** was

**€98.25**

**plus VAT &  
Shipping**



## Description

---

Pipette tips with filter, individually stuffed by machine, autoclavable. Particularly suitable for use in PCR technology and for work in microbiological and radioactive areas. Filter tips reduce the risk of cross-contamination and erroneous measurement results caused by contaminated pipette shafts and protect the shafts from contamination.

## Additional Information

---

No.	BN-2048
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Brand
old neoLab article no.	1-6546
EAN	4058072212834
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	200 µl
Volume MIN	5 µl
sterile	Yes
Pipette tips properties	pyrogen-free with certification with filter with box graduated RNase-free RNA-free PCR inhibitor-free endotoxin-free DNase-free DNA-free ATP-free
Fluid behavior	hydrophobic
Packing units tips	10 x 96
Model pipette tip	BRAND BIO-CERT Tips. Standard
Type certificate	Declaration of conformity (standardization)
Packaging type pipette tips (BRAND)	TipBox
for pipettes / dispensers of the manufacturers	Brand Gilson Thermo Scientific

