



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Filter Tips, ep Dualfilter T.I.P.S.® PCR clean, sterile, 0.2-5ml L, 175mm, purple

Special Price

€75.75 was

€101.00

**plus VAT &
Shipping**

Product Images



Description

ep Dualfilter T.I.P.S.®, PCR clean and sterile, 0.2 – 5 mL L, 175 mm, purple, colorless tips, 120 tips (5 racks × 24 tips). Replaces order no. 0030077725. – Two filter layers provide dual protection from aerosols¹⁾ and biomolecules to protect your pipette and sample – Free of PCR inhibitor additives for reliable results from downstream processes – Available in PCR clean/Sterile (sterile and pyrogen-free) – Constant control of each production batch by an independent laboratory guarantees highest purity and quality – Batch-specific certificates are available on the (URL: <https://www.eppendorf.com/certificates/>) [certificate homepage] – The certificate on filter efficacy according to EN 29463 is (URL: https://www.eppendorf.com/fileadmin/Main/05-Files/INT/Service_/_Support/CoQ_EPA_HEPA_Standard_of_epDualfilter_T.I.P.S._March_2020.pdf) [here] available – Modern, slim rack with resealable lid – Improved sustainability due to 19% lower polypropylene consumption per rack 1) An aerosol is a colloid of fine, solid particles or liquid droplets in air or another gas.

Additional Information

No.	VB-2638
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Eppendorf
EAN	4058072426187
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (violet)
Volume MAX	5 ml
Volume MIN	0.2 ml
sterile	Yes
Length	175 mm
Pipette tips properties	sterile PCR inhibitor-free
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	One rack/tray with box

