



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Filter Tips, ep Dualfilter T.I.P.S.® SealMax PCR clean, sterile, 2-100µl, 53mm, yellow

Special Price

**€144.75** was

**€193.00**

plus VAT &  
Shipping

### Product Images



## Description

ep Dualfilter T.I.P.S.® SealMax, PCR clean and sterile, 2 - 100 µL, 53 mm, yellow, colorless tips, 960 tips (10 racks × 96 tips). Replaces order no. 0030077822 - Reliable protection against liquids, aerosols<sup>1)</sup> and biomolecules - Certificate of filtration efficiency according to EN 29463 available at: (URL: <https://www.eppendorf.com/certificates>) [certificates homepage] - Free of PCR inhibitor additives - Protection of your pipette against aggressive liquids in case of overpipetting - Available in PCR clean/Sterile (Sterile and Pyrogen-free) - Constant control of each production batch by an independent laboratory ensures highest purity and quality - Batch specific certificates are available on (URL: <https://www.eppendorf.com/certificates/>) [certificates homepage] - Certificate on filter efficacy according to EN 29463 is (URL: [https://www.eppendorf.com/fileadmin/Main/05-Files/INT/Service\\_\\_\\_Support/CoQ\\_EPA\\_HEPA\\_Standard\\_of\\_epDualfilter\\_T.I.P.S.\\_March\\_2020.pdf](https://www.eppendorf.com/fileadmin/Main/05-Files/INT/Service___Support/CoQ_EPA_HEPA_Standard_of_epDualfilter_T.I.P.S._March_2020.pdf)) [here] available - Modern, slim rack with resealable lid - Improved sustainability due to 20% lower polypropylene consumption per rack 1) An aerosol is a colloid of fine, solid particles or liquid droplets in air or another gas.

## Additional Information

|                                       |                            |
|---------------------------------------|----------------------------|
| No.                                   | VB-2647                    |
| Typ Pipettenspitze                    | Filterspitzen              |
| Manufacturer (Brand)                  | Eppendorf                  |
| EAN                                   | 4058072426262              |
| Transport temperature                 | Room temperature           |
| autoclavable at 121°C                 | No                         |
| Color                                 | transparent (yellow)       |
| Volume MAX                            | 100 µl                     |
| Volume MIN                            | 2 µl                       |
| sterile                               | Yes                        |
| Length                                | 53 mm                      |
| Pipette tips properties               | sterile PCR inhibitor-free |
| Fluid behavior                        | hydrophobic                |
| Packaging type pipette tips (general) | One rack/tray with box     |

