



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf epDualfilter T.I.P.S. 1-10 ml L, individually sterile, 100 pcs./pack

Special Price

**€63.30** was

**€84.40**

**plus VAT &  
Shipping**

### Product Images



## Description

---

**Eppendorf epDualfilter T.I.P.S. are the first filter tips with contamination protection through a two-phase filter. The filter layers retain 100% of all aerosols and biomolecules.**

- double protection by two filter layers
- double protection for pipette and sample
- double protection against aerosols and biomolecules
- free of PCR inhibitor additives
- sterile, Eppendorf PCR clean and pyrogen-free
- IvD compliant
- batch certification

For use with Eppendorf Reference and Research as well as suitable for almost all common pipette models from other manufacturers.

**DNase-free, RNase-free, PCR inhibitor-free and human DNA-free.**

### **epT.I.P.S.<sup>®</sup> Racks**

Disposable boxes, not reusable, Order No. E-6521: 100 tips, individually packed

## Additional Information

No.	VB-0223
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6521
EAN	4058072405274
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (milky)
Volume MAX	10 ml
Volume MIN	1 ml
sterile	Yes
Pipette tips properties	with certification with filter with box graduated RNase-free PCR inhibitor-free DNase-free DNA-free
Fluid behavior	hydrophobic
Model pipette tip	Eppendorf ep Dualfilter T.I.P.S.
Type certificate	Batch certificate from the manufacturer
Packaging type Pipette tips (Eppendorf)	Racks
Packaging type pipette tips (general)	Carton. Tips individually packed
Packing units tips	100 x 1
for pipettes / dispensers of the manufacturers	Eppendorf

