



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf epDualfilter T.I.P.S. 100–5000 µl, sterile, in racks

Special Price

€76.50 was

€102.00

plus VAT &
Shipping

Product Images



Description

Eppendorf epDualfilter T.I.P.S. are the first filter tips with contamination protection through a two-phase filter. The filter layers retain 100% of all aerosols and biomolecules.

- double protection by two filter layers
- double protection for pipette and sample
- double protection against aerosols and biomolecules
- free of PCR inhibitor additives
- sterile, Eppendorf PCR clean and pyrogen-free
- IvD compliant
- batch certification

For use with Eppendorf Reference and Research as well as suitable for almost all common pipette models of other manufacturers.

DNase-free, RNase-free, PCR inhibitor-free and human DNA-free.

epT.I.P.S.[®] Racks

Disposable boxes, not reusable, Order No. E-6522: 5 boxes of 24 tips each

Additional Information

No.	VB-2637
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Eppendorf
old neoLab article no.	E-6522
EAN	4058072426170
Transport temperature	Room temperature
autoclavable at 121°C	No
Volume MAX	5 ml
Volume MIN	0.1 ml
sterile	Yes

