



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® InstaTip 200 µl filter pipette tips,  
sterile, DNase and RNase free, 10 x 96**

**€138.60  
plus VAT &  
Shipping**

**Product Images**

---

## Description

---

InstaTip 1000 µl Filter PipetteTip, sterilized by irradiation, DNase & RNase free

## Additional Information

---

No.	7-2005
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	neoLab
EAN	4066292687646
Transport temperature	Room temperature
Color	transparent (clear)
Volume MAX	200 µl
Volume MIN	1 µl
sterile	Yes
Pipette tips properties	sterile nose-free Rnose-free
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	Multiple racks/ trays each with box

