



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® 10µl Filter Tips, Sterile, DNase and RNase free, PP

€49.90
plus VAT &
Shipping

Product Images



Description

neolab filter tips offer you optimum protection against aerosols and splashes for your samples and pipettes. The filter tips are particularly suitable for liquid handling in the PCR technique. In addition, they are used in all areas where either the sample or the pipette needs to be protected during pipetting (e.g. microbiology, working with radioactive substances, working with blood).

The hydrophobic filter prevents the transfer of liquids and aerosols between the sample and the shaft of the pipette used. This also reduces the risk of cross-contamination during work.

The **additional extra-long version** allows safe and clean pipetting with narrow vessels. The filter tips are compatible with all common models, e.g. pipettes from Eppendorf, Brand and Gilson.

- DNase- and RNase-free
- sterile
- conveniently stored in box with color-coded trays

Additional Information

No.	9-2000
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	neolabLine
EAN	4058072983765
Transport temperature	Room temperature
diameter	5,8 mm
Color	transparent (clear)
Volume MAX	10 µl
Volume MIN	0.1 µl
sterile	Yes
Length	32 mm
Pipette tips properties	sterile nose-free RNase-free
Color code	white
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	One rack/tray with box

