



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® 200µl Filter Tips, Sterile, DNase and RNase free, PP

**€46.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**neolab filter tips** offer you optimum protection against aerosols and splashes for your samples and pipettes. The filter tips are particularly suitable for liquid handling in the PCR technique. In addition, they are used in all areas where either the sample or the pipette needs to be protected during pipetting (e.g. microbiology, working with radioactive substances, working with blood).

The hydrophobic filter prevents the transfer of liquids and aerosols between the sample and the shaft of the pipette used. This also reduces the risk of cross-contamination during work.

The **additional extra-long version** allows safe and clean pipetting with narrow vessels. The filter tips are compatible with all common models, e.g. pipettes from Eppendorf, Brand and Gilson.

- DNase- and RNase-free
- sterile
- conveniently stored in box with color-coded trays

## Additional Information

---

|                                       |                        |
|---------------------------------------|------------------------|
| No.                                   | 9-2003                 |
| Typ Pipettenspitze                    | Filterspitzen          |
| Manufacturer (Brand)                  | neolabLine             |
| EAN                                   | 4058072983796          |
| Transport temperature                 | Room temperature       |
| diameter                              | 7,5 mm                 |
| Color                                 | transparent (clear)    |
| Volume MAX                            | 200 µl                 |
| Volume MIN                            | 10 µl                  |
| sterile                               | Yes                    |
| Length                                | 50.5 mm                |
| Color code                            | bordeaux red           |
| Fluid behavior                        | hydrophobic            |
| Packaging type pipette tips (general) | One rack/tray with box |

