



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® 10µl Filter Tips, Low retention, Sterile, DNase and RNase free, PP, Extra-long

**€30.50**  
**plus VAT &**  
**Shipping**

### Product Images



# Description

**neolab filter tips** offer you optimum protection from aerosols and splashes for your samples and pipettes. The filter tips are particularly suitable for liquid handling in the PCR technique. In addition, they are used in all areas where either the sample or the pipette needs to be protected during pipetting (e.g. microbiology, working with radioactive substances, working with blood).

The hydrophobic filter prevents the transfer of liquids and aerosols between the sample and the shaft of the pipette used. This also reduces the risk of cross-contamination during work.

The **additional extra-long version** allows safe and clean pipetting with narrow vessels. The filter tips are compatible with all common models, e.g. pipettes from Eppendorf, Brand and Gilson.

- DNase- and RNase-free
- sterile
- conveniently stored in box with color-coded trays

# Additional Information

No.	9-2010
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	neolabLine
EAN	4058072983864
Transport temperature	Room temperature
diameter	5,8 mm
Color	transparent (clear)
Volume MAX	10 µl
Volume MIN	0.1 µl
sterile	Yes
Length	45.8 mm
Fluid behavior	hydrophobic
Packaging type pipette tips (general)	One rack/tray with box

