



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Gilson Blister Refill DFL10ST, 0.1-10 µl, sterile, 10 x 96 pcs/pack

**€132.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**The refill variant for Gilson PIPETMAN® Filtertips Diamond®, sterile, in blister packaging.**

Filtertips can help prevent aerosol contamination.

Gamma-sterilized, free of detectable RNases, DNases, DNA, RNA, and ATP.

Individually packaged racks with color coding compatible with Gilson pipette models.

For maximum cleanliness and ease of use.

Simply remove the cover foil from the packaging and transfer the ready-to-use rack with tips into the corresponding box.

**Blister Refill™**

Sterile, color-coded sealed.

10 trays individually packed with 96 tips each

## Additional Information

No.	CF-0445
Typ Pipettenspitze	Filterspitzen
Manufacturer (Brand)	Gilson
old neoLab article no.	7-4464
EAN	4058072180492
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Volume MAX	10 µl
Volume MIN	0.1 µl
sterile	Yes
Pipette tips properties	pyrogen-free with certification with filter with box extra long RNase-free RNA-free DNase-free DNA-free ATP-free
Fluid behavior	hydrophobic
Designation tip (Gilson)	DFL10ST
Packaging type pipette tips (general)	Multiple racks/ trays each with box
Packing units tips	10 x 96
Model pipette tip	Gilson PIPETMAN Tips Diamond
Type certificate	Sterilization certificate from the manufacturer
Packaging type pipette tips (Gilson)	Blister
for pipettes / dispensers of the manufacturers	Gilson

